

Patient Prep Bag Content:

2 each of 2-inch medical tape, Bio-Hazard bag, 3 pairs of large gloves, 4 of triangular bandages.



Tie the backboard patient kit to the backboard as shown here.



Backboard Setup

Backboards need to be checked for contents daily.

Secure the backboard bag to the backboard.



Place a piece of tape indicating "complete", date, and the person checking. Return the backboard to the wall storage



Patient Prep Bag



2 each of head-packs, 1 each of Spider Strap, 1 adult c-collar, 1 pediatric c-collar. 3-4 pcs of carpet padding, and one patient prep bag.